

ALTAN-Medical-Service - Zweifelsheimer Weg 19/b - 91448 Emskirchen

ALTAN-Medical-Service



Zweifelsheimer Weg 19/b
91448 Emskirchen
Tel.: +49 (0) 9101-9055530
Fax.: +49 (0) 9101-9055539
E-Mail: info@altan-medical.com
Internet: www.altan-medical.com
VAT-No. DE273938430

Sehr geehrte Damen und Herren,
Dear Sir and Madam,

wir können Ihnen die angefragte Maschine - Zwischenverkauf/Irrtümer vorbehalten - unverbindlich wie folgt anbieten:
We can offer you the following requested machine(s) - subject to prior sale:

Anesthesia - Anesthesiology

Anästhesiegerät - Beatmungsgerät / Anesthesiology Equipment DRÄGER ZEUS SW 3.08

Year 2008 - Serial-No. ARZF-0020

Lung - protective ventilation applications

2 DIVA Target Controlled Anaesthesia - dosing module Sevoflurane - Desflurane

Drägersorb CLIC Drägersorb is high-quality soda lime developed by one of the leading manufacturers of anaesthesia equipment

Dräger SmartPilot® View with Touchscreen DRÄGER Infinity C500

Dräger VentStar® Anaesthesia Breathing Circuits

Technical Data:

Gas delivery

Modes of operation Fresh-gas control, auto control with insp. O2 control

Fresh-gas flow 0.25 ? 18 l/min (fresh-gas control); 0 ? 18 l/min (auto control); closed system mode (uptake)

O2 control 21 ? 100 % (in fresh-gas control with AIR); Oxygen Ratio Controller

Carrier gases N2O and Air

O2 flush > 35 l/min

O2 safety flow 0 ? 12 l/min

Auxiliary O2 flow tube 0 ? 16 l/min (e.g. for regional anaesthesia)

Anaesthetic agent delivery DIVA (Direct Injection of Volatile Agents)

Mode of operation of DIVA modules Fresh-gas control, auto control with expiratory agent control (2 plug-in slots for ISO, SEV, DES)

Max. filling volume 315 ml

Delivery range: Fresh-gas control (fresh-gas settings) ISO: 0 ? 5 kPa; SEV: 0 ? 8 kPa; DES: 0 ? 18 kPa

Delivery range: Auto control (expiratory agent settings) ISO: 0 ? 2.5 kPa; SEV: 0 ? 5 kPa;

You are welcome to contact us by Messenger (WhatsApp-Telegram-Viber-Signal): +49(0)178-3648464

Visit our website for the latest used machinery: <https://altan-medical.com/machines-list>

DES: 0 ? 12 kPa

General Technical Data:

Data exchange up to 6 USB ports

Serial interface 2 x RS232

Ethernet interface (1) For access to hospital network, network printer, or remote inspection

Ethernet interface (2) For access to hospital network printer, or remote inspection

Ethernet interface (3) Infinity network connection

Power supply 100 ? 240 V, 45 ? 65 Hz

Battery backup min. 30, typically 90 minutes

Dimensions (H x W x D) 165 x 110 x 70 cm

Weight 185 kg incl. 2 DIVA module

Online-Video-Inspection by Skype-Video

We would be very pleased with your visit - more machines on Stock

Available Immediately - Can be inspect

Standort / Location: On Stock Emskirchen / Nürnberg

Price: On request Euro FCA Emskirchen/ Nürnberg

[View: Machine on Website](#)

You are welcome to contact us by Messenger (WhatsApp-Telegram-Viber-Signal): +49(0)178-3648464

Visit our website for the latest used machinery: <https://altan-medical.com/machines-list>

ALTAN-Medical-Service
Zweifelsheimer Weg 19/b
91448 Emskirchen
Tel.: +49 (0) 9101-9055530

Bankverbindung
Commerzbank Nürnberg
IBAN: DE10 7608 0040 0152 8774 00
BIC: DRES DE FF760


Geschäftsführer Altan
Amtsgericht Fürth - HRB-No. 17959
Steuer-Nr. 203/126/00529
Es gelten unsere AGB.

ALTAN-Medical-Service - Zweifelsheimer Weg 19/b - 91448 Emskirchen



ALTAN-Medical-Service

Zweifelsheimer Weg 19/b
91448 Emskirchen
Tel.: +49 (0) 9101-9055530
Fax.: +49 (0) 9101-9055539
E-Mail: [info\(at\)altan-medical.com](mailto:info(at)altan-medical.com)
Internet: www.altan-medical.com
VAT-No. DE273938430



You are welcome to contact us by Messenger (WhatsApp-Telegram-Viber-Signal): +49(0)178-3648464
Visit our website for the latest used machinery: <https://altan-medical.com/machines-list>